

# GOWNS



## Surgical Gown (standard performance)

### • Product description


- Made of Blue alcohol-repellent SMMS Non-woven fabric , it was designed to meet standard performance for lower-fluid.

### • Product features

- alcohol-repellent anti-bacteria.
- Safe and comfortable,
- Comply with EN13795:standard performance.
- Comply with AAMIPB70:level 2 level 3

### • Application

It is intended to be worn by operation room personnel during surgical procedures to protect the surgical patient and the operating room personnel from the transfer of microorganisms, body fluids.



### • Package

Size:	M / L / LL / XL / XLL / XXL
Colors:	Blue
Pack:	